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- Eguchi Y et al: Prevalence and associated metabolic factors of nonalcoholic fatty liver disease in the general population from 2009 to 2010 in Japan: a multicenter large retrospective study. J Gastroenterol 47 (5): 586-95, 2012
- Joseph AEA et al: Ultrasound in the detection of chronic liver disease (the "bright liver"). Brit J Radiol 52 (615): 184-8, 1979
- Yajima Y et al: Ultrasonographical diagnosis of fatty liver: Significance of the liver -kidney contrast. Tohoku J exp Med 139 (1): 43-50, 1983
- Kleiner DE et al: Design and validation of a histological scoring system for nonalcoholic fatty liver disease. Hepatology 41 (6): 1313-21, 2005
- Targher G et al: Nonalcoholic fatty liver disease is independently associated with an increased incidence of cardiovascular events in type 2 diabetic patients. Diabetes Care 30 (8): 2119-21, 2007
- Urena M et al: An approach to the rational use of steatotic donor livers in liver transplantation. Hepato-Gastroenterol 46 (26): 1164-73, 1999
- Kudo M et al: JSUM ultrasound elastography practice guidelines: liver. J Med Ultrason (2001) 40 (4): 325-57, 2013
- Herrmann E et al: Assessment of biopsy-proven liver fibrosis by two-dimensional shear wave elastography: An individual patient data-based meta-analysis. Hepatology 67 (1): 260-272, 2018
- Angulo P et al: Liver Fibrosis, but No Other Histologic Features, Is Associated With Long-term Outcomes of Patients With Nonalcoholic Fatty Liver Disease. Gastroenterology 149 (2): 389-97, 2015
- Tada T et al: Progression of liver fibrosis is associated with non-liver-related mortality in patients with nonalcoholic fatty liver disease. Hepatol Commun 1 (9): 899-910, 2017
- European Association for Study of Liver et al: EASL-ALEH Clinical Practice Guidelines: Non-invasive tests for evaluation of liver disease severity and prognosis. J Hepatol 63 (1): 237-64, 2015

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